

# RAB1B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00675-1

## Capture Antibody Information

<b>Catalog Number:</b> 83687-2-PBS	<b>Clone ID:</b> 240635G11	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> RAB1B, member RAS oncogene family
<b>Isotype:</b> IgG	<b>GenBank:</b> BC071169	<b>Gene ID:</b> 81876
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag12197	

## Detection Antibody Information

<b>Catalog Number:</b> 83687-1-PBS	<b>Clone ID:</b> 240635B3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> RAB1B, member RAS oncogene family
<b>Isotype:</b> IgG	<b>GenBank:</b> BC071169	<b>Gene ID:</b> 81876
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag12197	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.781-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP00675-1 targets RAB1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RAB1B Recombinant antibody, PBS Only (Capture) 83687-2-PBS (240635G11). 100 µg. Concentration 1 mg/mL.

Detection antibody: RAB1B Recombinant antibody, PBS Only (Capture/Detector) 83687-1-PBS (240635B3). 100 µg. Concentration 1 mg/mL.

Alternative RAB1B matched antibody pairs: MP00675-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00675-1, RAB1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83687-2-PBS. Detection antibody: 83687-1-PBS. Standard: Ag12197. Range: 0.781-100 ng/mL.