

IL-7Ra/CD127 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00704-1

Capture Antibody Information

Catalog Number: 83723-2-PBS	Clone ID: 240749C11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: interleukin 7 receptor
Isotype: IgG	GenBank: AAC83204.1	Gene ID: 3575
Purification Method: Protein A purification		

Detection Antibody Information

Catalog Number: 83723-1-PBS	Clone ID: 240749B2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: interleukin 7 receptor
Isotype: IgG	GenBank: AAC83204.1	Gene ID: 3575
Purification Method: Protein A purification		

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00704-1 targets IL-7Ra/CD127 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL-7Ra/CD127 Recombinant antibody, PBS Only (Capture/Detector) 83723-2-PBS (240749C11). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-7Ra/CD127 Recombinant antibody, PBS Only (Detector) 83723-1-PBS (240749B2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00704-1, IL-7Ra/CD127 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83723-2-PBS. Detection antibody: 83723-1-PBS. Standard: Eg1279. Range: 0.625-80 ng/mL