

# ACOX1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00709-1

## Capture Antibody Information

**Catalog Number:**  
83731-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240754G9  
**Reactivity:**  
human  
**GenBank:**  
BC010425  
**Immunogen Catalog Number:**  
Ag1405

**Conjugate:**  
Unconjugated  
**Full name:**  
acyl-Coenzyme A oxidase 1, palmitoyl  
**Gene ID:**  
51

## Detection Antibody Information

**Catalog Number:**  
83731-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240754D4  
**Reactivity:**  
human  
**GenBank:**  
BC010425  
**Immunogen Catalog Number:**  
Ag1405

**Conjugate:**  
Unconjugated  
**Full name:**  
acyl-Coenzyme A oxidase 1, palmitoyl  
**Gene ID:**  
51

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.78-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00709-1 targets ACOX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AOX Recombinant antibody, PBS Only (Capture) 83731-4-PBS (240754G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: AOX Recombinant antibody, PBS Only (Capture/Detector) 83731-3-PBS (240754D4). 100 µg. Concentration 1 mg/mL.

Alternative ACOX1 matched antibody pairs: MP00709-2, MP00709-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00709-1, AOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83731-4-PBS. Detection antibody: 83731-3-PBS. Standard: Ag1405. Range: 0.78-100 ng/mL