

VEGF-C Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00710-2**

Capture Antibody Information

Catalog Number:

83727-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240772C1

Reactivity:

human

GenBank:

BC035212

Conjugate:

Unconjugated

Full name:

vascular endothelial growth factor C

Gene ID:

7424

Detection Antibody Information

Catalog Number:

83727-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240772A5

Reactivity:

human

GenBank:

BC035212

Conjugate:

Unconjugated

Full name:

vascular endothelial growth factor C

Gene ID:

7424

Applications

Tested Applications:

Cytometric bead array

Range:

0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00710-2 targets VEGF-C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VEGF-C Recombinant antibody, PBS Only (Capture/Detector) 83727-3-PBS (240772C1). 100 µg. Concentration 1 mg/mL.

Detection antibody: VEGF-C Recombinant antibody, PBS Only (Detector) 83727-1-PBS (240772A5). 100 µg. Concentration 1 mg/mL.

Alternative VEGF-C matched antibody pairs: MP00710-1, MP00710-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00710-2, VEGF-C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83727-3-PBS. Detection antibody: 83727-1-PBS. Standard: Eg32025. Range: 0.313-20 ng/mL