

Mouse Apoe Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00712-2

Capture Antibody Information

Catalog Number:	Clone ID:	Conjugate:
83728-2-PBS	240786G1	Unconjugated
Host:	Reactivity:	Full name:
Rabbit	mouse	Apolipoprotein E
Isotype:	GenBank:	Gene ID:
IgG	NM_001305819.1	11816
Purification Method:	Immunogen Catalog Number:	
Protein A purification	Eg1093	

Detection Antibody Information

Catalog Number:	Clone ID:	Conjugate:
83728-1-PBS	240786C6	Unconjugated
Host:	Reactivity:	Full name:
Rabbit	mouse	Apolipoprotein E
Isotype:	GenBank:	Gene ID:
IgG	NM_001305819.1	11816
Purification Method:	Immunogen Catalog Number:	
Protein A purification	Eg1093	

Applications

Tested Applications:	Range:	Recommended Dilutions:
Cytometric bead array	0.313-40 ng/mL (Cytometric Bead Array)	It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00712-2 targets Apoe in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Apolipoprotein E Recombinant antibody, PBS Only 83728-2-PBS (240786G1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Mouse Apolipoprotein E Recombinant antibody, PBS Only (Detector) 83728-1-PBS (240786C6). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00712-2, MOUSE Apoe Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83728-2-PBS. Detection antibody: 83728-1-PBS. Standard: Eg1093. Range: 0.313-40 ng/mL