

# CBR1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00724-1

## Capture Antibody Information

Catalog Number: 83736-2-PBS	Clone ID: 240707F4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: carbonyl reductase 1
Isotype: IgG	GenBank: BC002511	Gene ID: 873
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34830	

## Detection Antibody Information

Catalog Number: 83736-1-PBS	Clone ID: 240707A7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: carbonyl reductase 1
Isotype: IgG	GenBank: BC002511	Gene ID: 873
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34830	

## Applications

Tested Applications: Cytometric bead array	Range: 1.56-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

## Product Information

MP00724-1 targets CBR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CBR1 Recombinant antibody, PBS Only (Capture) 83736-2-PBS (240707F4). 100 µg. Concentration 1 mg/mL.

Detection antibody: CBR1 Recombinant antibody, PBS Only (Detector) 83736-1-PBS (240707A7). 100 µg. Concentration 1 mg/mL.

Alternative CBR1 matched antibody pairs: MP00724-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00724-1, CBR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83736-2-PBS. Detection antibody: 83736-1-PBS. Standard: Ag34830. Range: 1.56-100 ng/mL