

YAF2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00726-2**

Capture Antibody Information

Catalog Number: 83798-4-PBS	Clone ID: 240834B12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: YY1 associated factor 2
Isotype: IgG	GenBank: BC037777	Gene ID: 10138
Purification Method: Protein A purification	Immunogen Catalog Number: Ag5023	

Detection Antibody Information

Catalog Number: 83798-2-PBS	Clone ID: 240834F11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: YY1 associated factor 2
Isotype: IgG	GenBank: BC037777	Gene ID: 10138
Purification Method: Protein A purification	Immunogen Catalog Number: Ag5023	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00726-2 targets YAF2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: YAF2 Recombinant antibody, PBS Only (Capture) 83798-4-PBS (240834B12). 100 µg. Concentration 1 mg/mL.

Detection antibody: YAF2 Recombinant antibody, PBS Only (Detector) 83798-2-PBS (240834F11). 100 µg. Concentration 1 mg/mL.

Alternative YAF2 matched antibody pairs: MP00726-1, MP00726-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00726-2, YAF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83798-4-PBS. Detection antibody: 83798-2-PBS. Standard: Ag5023. Range: 0.625-40 ng/mL