

# CTLA-4/CD152 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00736-2**

## Capture Antibody Information

**Catalog Number:**  
83778-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240793F6  
**Reactivity:**  
human  
**GenBank:**  
BC070162

**Conjugate:**  
Unconjugated  
**Full name:**  
cytotoxic T-lymphocyte-associated protein 4  
**Gene ID:**  
1493

## Detection Antibody Information

**Catalog Number:**  
83778-5-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240793F7  
**Reactivity:**  
human  
**GenBank:**  
BC070162

**Conjugate:**  
Unconjugated  
**Full name:**  
cytotoxic T-lymphocyte-associated protein 4  
**Gene ID:**  
1493

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00736-2 targets CTLA-4/CD152 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CTLA-4/CD152 Recombinant antibody, PBS Only (Capture) 83778-4-PBS (240793F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: CTLA-4/CD152 Recombinant antibody, PBS Only (Detector) 83778-5-PBS (240793F7). 100 µg. Concentration 1 mg/mL.

Alternative CTLA-4/CD152 matched antibody pairs: MP00736-1, MP00736-3, MP00736-4, MP00907-1, MP00907-2, MP00907-3, MP00907-4, MP50106-1, MP50106-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00736-2, CTLA4/CD152 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83778-4-PBS. Detection antibody: 83778-5-PBS. Standard: Eg0972. Range: 0.313-40 ng/mL.