

BCL9 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00739-1

Capture Antibody Information

Catalog Number: 83760-1-PBS	Clone ID: 240682B4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: B-cell CLL/lymphoma 9
Isotype: IgG	GenBank: BC116451	Gene ID: 607
Purification Method: Protein A purification	Immunogen Catalog Number: Ag19158	

Detection Antibody Information

Catalog Number: 83760-2-PBS	Clone ID: 240682G8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: B-cell CLL/lymphoma 9
Isotype: IgG	GenBank: BC116451	Gene ID: 607
Purification Method: Protein A purification	Immunogen Catalog Number: Ag19158	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
-----------------------------------------------	--------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP00739-1 targets BCL9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCL9 Recombinant antibody, PBS Only (Capture) 83760-1-PBS (240682B4). 100 μ g. Concentration 1 mg/mL.

Detection antibody: BCL9 Recombinant antibody, PBS Only (Detector) 83760-2-PBS (240682G8). 100 μ g. Concentration 1 mg/mL.

Alternative BCL9 matched antibody pairs: MP00739-2, MP00739-3, MP00739-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00739-1, BCL9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83760-1-PBS. Detection antibody: 83760-2-PBS. Standard: Ag19158. Range: 0.313-40 ng/mL