

MAFB Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00750-1

Capture Antibody Information

Catalog Number:
83783-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240821D5

Reactivity:
Human

GenBank:
BC036689

Immunogen Catalog Number:
Ag14102

Conjugate:
Unconjugated

Full name:
v-maf musculoaponeurotic
fibrosarcoma oncogene homolog B
(avian)

Gene ID:
9935

Detection Antibody Information

Catalog Number:
83783-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240821F12

Reactivity:
Human

GenBank:
BC036689

Immunogen Catalog Number:
Ag14102

Conjugate:
Unconjugated

Full name:
v-maf musculoaponeurotic
fibrosarcoma oncogene homolog B
(avian)

Gene ID:
9935

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00750-1 targets MAFB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MAFB Recombinant antibody, PBS Only (Capture) 83783-2-PBS (240821D5). 100 µg. Concentration 1 mg/mL.

Detection antibody: MAFB Recombinant antibody, PBS Only (Detector) 83783-4-PBS (240821F12). 100 µg. Concentration 1 mg/mL.

Alternative MAFB matched antibody pairs: MP00750-2, MP00750-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00750-1, MAFB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83783-2-PBS. Detection antibody: 83783-4-PBS. Standard: Ag14102. Range: 0.313-40 ng/mL