

TNC/Tenascin-C Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00755-1

Capture Antibody Information

Catalog Number: 83767-2-PBS	Clone ID: 240844D8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: tenascin C
Isotype: IgG	GenBank: BC151843	Gene ID: 3371
Purification Method: Protein A purification	Immunogen Catalog Number: Ag27122	

Detection Antibody Information

Catalog Number: 83767-4-PBS	Clone ID: 240844F12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: tenascin C
Isotype: IgG	GenBank: BC151843	Gene ID: 3371
Purification Method: Protein A purification	Immunogen Catalog Number: Ag27122	

Applications

Tested Applications: Cytometric bead array	Range: 0.78-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

Product Information

MP00755-1 targets TNC/Tenascin-C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNC/Tenascin-C Recombinant antibody, PBS Only (Capture) 83767-2-PBS (240844D8). 100 µg. Concentration 1 mg/mL.

Detection antibody: TNC/Tenascin-C Recombinant antibody, PBS Only (Detector) 83767-4-PBS (240844F12). 100 µg. Concentration 1 mg/mL.

Alternative TNC/Tenascin-C matched antibody pairs: MP00755-2, MP00755-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00755-1, TNC/Tenascin-C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83767-2-PBS. Detection antibody: 83767-4-PBS. Standard: Ag27122. Range: 0.78-100 ng/mL.