

HNRNPA2B1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00768-2

Capture Antibody Information

Catalog Number: 83773-1-PBS	Clone ID: 240831E9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: heterogeneous nuclear ribonucleoprotein A2/B1
Isotype: IgG	GenBank: BC000506	Gene ID: 3181
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6563	

Detection Antibody Information

Catalog Number: 83773-3-PBS	Clone ID: 240831A8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: heterogeneous nuclear ribonucleoprotein A2/B1
Isotype: IgG	GenBank: BC000506	Gene ID: 3181
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6563	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

Product Information

MP00768-2 targets HNRNPA2B1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HNRNPA2B1 Recombinant antibody, PBS Only (Capture) 83773-1-PBS (240831E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: HNRNPA2B1 Recombinant antibody, PBS Only (Detector) 83773-3-PBS (240831A8). 100 µg. Concentration 1 mg/mL.

Alternative HNRNPA2B1 matched antibody pairs: MP00768-1, MP00768-3, MP00768-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00768-2, HNRNPA2B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83773-1-PBS. Detection antibody: 83773-3-PBS. Standard: Ag6563. Range: 0.313-40 ng/mL