

BRUNOL5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00774-1

Capture Antibody Information

Catalog Number:
83789-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240693E4

Reactivity:
Human

GenBank:
BC028101

Immunogen Catalog Number:
Ag4416

Conjugate:
Unconjugated

Full name:
bruno-like 5, RNA binding protein (Drosophila)

Gene ID:
60680

Detection Antibody Information

Catalog Number:
83789-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240693G5

Reactivity:
Human

GenBank:
BC028101

Immunogen Catalog Number:
Ag4416

Conjugate:
Unconjugated

Full name:
bruno-like 5, RNA binding protein (Drosophila)

Gene ID:
60680

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00774-1 targets BRUNOL5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BRUNOL5 Recombinant antibody, PBS Only (Capture) 83789-3-PBS (240693E4). 100 µg. Concentration 1 mg/mL.

Detection antibody: BRUNOL5 Recombinant antibody, PBS Only (Detector) 83789-2-PBS (240693G5). 100 µg. Concentration 1 mg/mL.

Alternative BRUNOL5 matched antibody pairs: MP00774-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00774-1, BRUNOL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83789-3-PBS. Detection antibody: 83789-2-PBS. Standard: Ag4416. Range: 0.625-80 ng/mL