

# HIF2 $\alpha$ /EPAS1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00775-1

## Capture Antibody Information

**Catalog Number:**  
83790-4-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240736H7

**Reactivity:**  
human

**GenBank:**  
BC051338

**Immunogen Catalog Number:**  
Ag15199

**Conjugate:**  
Unconjugated

**Full name:**  
endothelial PAS domain protein 1

**Gene ID:**  
2034

## Detection Antibody Information

**Catalog Number:**  
83790-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240736F2

**Reactivity:**  
human, rat

**GenBank:**  
BC051338

**Immunogen Catalog Number:**  
Ag15199

**Conjugate:**  
Unconjugated

**Full name:**  
endothelial PAS domain protein 1

**Gene ID:**  
2034

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.78-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00775-1 targets HIF2  $\alpha$  /EPAS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HIF2  $\alpha$  /EPAS1 Recombinant antibody, PBS Only (Capture) 83790-4-PBS (240736H7). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: HIF2  $\alpha$  /EPAS1 Recombinant antibody, PBS Only (Detector) 83790-1-PBS (240736F2). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00775-1, HIF2  $\alpha$  /EPAS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83790-4-PBS. Detection antibody: 83790-1-PBS. Standard: Ag15199. Range: 0.78-100 ng/mL.