

HIF2 α /EPAS1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00775-3

Capture Antibody Information

Catalog Number:
83790-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240736H7

Reactivity:
human

GenBank:
BC051338

Immunogen Catalog Number:
Ag15199

Conjugate:
Unconjugated

Full name:
endothelial PAS domain protein 1

Gene ID:
2034

Detection Antibody Information

Catalog Number:
83790-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240736F12

Reactivity:
human

GenBank:
BC051338

Immunogen Catalog Number:
Ag15199

Conjugate:
Unconjugated

Full name:
endothelial PAS domain protein 1

Gene ID:
2034

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00775-3 targets HIF2 α /EPAS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HIF2 α /EPAS1 Recombinant antibody, PBS Only (Capture) 83790-4-PBS (240736H7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: HIF2 α /EPAS1 Recombinant antibody, PBS Only (Detector) 83790-3-PBS (240736F12). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00775-3, HIF2 α /EPAS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83790-4-PBS. Detection antibody: 83790-3-PBS. Standard: Ag15199. Range: 1.56-100 ng/mL.