

EXTL2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00784-1

Capture Antibody Information

| | | |
|--|-------------------------------------|---|
| Catalog Number: 83821-2-PBS | Clone ID: 240919G2 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: exostoses (multiple)-like 2 |
| Isotype: IgG | GenBank: BC036015 | Gene ID: 2135 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag3548 | |

Detection Antibody Information

| | | |
|--|-------------------------------------|---|
| Catalog Number: 83821-1-PBS | Clone ID: 240919B10 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: exostoses (multiple)-like 2 |
| Isotype: IgG | GenBank: BC036015 | Gene ID: 2135 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag3548 | |

Applications

| | | |
|---|--|--|
| Tested Applications: Cytometric bead array | Range: 0.313-20 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|---|--|--|

Product Information

MP00784-1 targets EXTL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EXTL2 Recombinant antibody, PBS Only (Capture) 83821-2-PBS (240919G2). 100 µg. Concentration 1 mg/mL.

Detection antibody: EXTL2 Recombinant antibody, PBS Only (Detector) 83821-1-PBS (240919B10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00784-1, EXTL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83821-2-PBS. Detection antibody: 83821-1-PBS. Standard: Ag3548. Range: 0.313-20 ng/mL