

MIOS Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00801-1

Capture Antibody Information

Catalog Number:
83856-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240832A10

Reactivity:
human

GenBank:
BC005883

Immunogen Catalog Number:
Ag15148

Conjugate:
Unconjugated

Full name:
missing oocyte, meiosis regulator, homolog (Drosophila)

Gene ID:
54468

Detection Antibody Information

Catalog Number:
83856-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240832D7

Reactivity:
human

GenBank:
BC005883

Immunogen Catalog Number:
Ag15148

Conjugate:
Unconjugated

Full name:
missing oocyte, meiosis regulator, homolog (Drosophila)

Gene ID:
54468

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00801-1 targets MIOS in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MIOS Recombinant antibody, PBS Only (Capture/Detector) 83856-1-PBS (240832A10). 100 µg. Concentration 1 mg/mL.

Detection antibody: MIOS Recombinant antibody, PBS Only (Detector) 83856-2-PBS (240832D7). 100 µg. Concentration 1 mg/mL.

Alternative MIOS matched antibody pairs: MP00801-2, MP00801-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00801-1, MIO5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83856-1-PBS. Detection antibody: 83856-2-PBS. Standard: Ag15148. Range: 0.625-80 ng/mL