

CNNM4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00802-1

Capture Antibody Information

Catalog Number:
83824-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240873B3

Reactivity:
human

GenBank:
BC063295

Immunogen Catalog Number:
Ag5121

Conjugate:
Unconjugated

Full name:
cyclin M4

Gene ID:
26504

Detection Antibody Information

Catalog Number:
83824-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240873D3

Reactivity:
human

GenBank:
BC063295

Immunogen Catalog Number:
Ag5121

Conjugate:
Unconjugated

Full name:
cyclin M4

Gene ID:
26504

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00802-1 targets CNNM4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CNNM4 Recombinant antibody, PBS Only (Capture) 83824-1-PBS (240873B3). 100 µg. Concentration 1 mg/mL.

Detection antibody: CNNM4 Recombinant antibody, PBS Only (Detector) 83824-2-PBS (240873D3). 100 µg. Concentration 1 mg/mL.

Alternative CNNM4 matched antibody pairs: MP00802-2, MP00802-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00802-1, CNNM4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83824-1-PBS. Detection antibody: 83824-2-PBS. Standard: Ag5121. Range: 0.313-40 ng/mL