

SLC9A3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00818-2

Capture Antibody Information

Catalog Number:
83832-5-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240999F2

Reactivity:
human

GenBank:
BC101669

Immunogen Catalog Number:
Ag25637

Conjugate:
Unconjugated

Full name:
solute carrier family 9 (sodium/hydrogen exchanger), member 3

Gene ID:
6550

Detection Antibody Information

Catalog Number:
83832-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240999B6

Reactivity:
human

GenBank:
BC101669

Immunogen Catalog Number:
Ag25637

Conjugate:
Unconjugated

Full name:
solute carrier family 9 (sodium/hydrogen exchanger), member 3

Gene ID:
6550

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00818-2 targets SLC9A3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NHE3 Recombinant antibody, PBS Only (Capture) 83832-5-PBS (240999F2). 100 µg. Concentration 1 mg/mL.

Detection antibody: NHE3 Recombinant antibody, PBS Only (Detector) 83832-3-PBS (240999B6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00818-2, SLC9A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83832-5-PBS. Detection antibody: 83832-3-PBS. Standard: Ag25637. Range: 0.625-40 ng/mL