

USP53 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00821-3**

Capture Antibody Information

Catalog Number:
83846-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240971E1
Reactivity:
human
GenBank:
NM_019050
Immunogen Catalog Number:
Ag28923

Conjugate:
Unconjugated
Full name:
ubiquitin specific peptidase 53
Gene ID:
54532

Detection Antibody Information

Catalog Number:
83846-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240971H8
Reactivity:
human
GenBank:
NM_019050
Immunogen Catalog Number:
Ag28923

Conjugate:
Unconjugated
Full name:
ubiquitin specific peptidase 53
Gene ID:
54532

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00821-3 targets USP53 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: USP53 Recombinant antibody, PBS Only (Capture/Detector) 83846-2-PBS (240971E1). 100 µg. Concentration 1 mg/mL.

Detection antibody: USP53 Recombinant antibody, PBS Only (Detector) 83846-3-PBS (240971H8). 100 µg. Concentration 1 mg/mL.

Alternative USP53 matched antibody pairs: MP00821-1, MP00821-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00821-3, USP53 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83846-2-PBS. Detection antibody: 83846-3-PBS. Standard: Ag28923. Range: 1.56-100 ng/mL