

CXCL8/IL-8 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00830-2

Capture Antibody Information

Catalog Number:
83903-2-PBSHost:
RabbitIsotype:
IgGPurification Method:
Protein A purificationClone ID:
240889F6Reactivity:
humanGenBank:
BC013615Conjugate:
UnconjugatedFull name:
interleukin 8Gene ID:
3576

Detection Antibody Information

Catalog Number:
83903-1-PBSHost:
RabbitIsotype:
IgGPurification Method:
Protein A purificationClone ID:
240889C2Reactivity:
humanGenBank:
BC013615Conjugate:
UnconjugatedFull name:
interleukin 8Gene ID:
3576

Applications

Tested Applications:
Cytometric bead arrayRange:
0.625-40 ng/mL (Cytometric Bead Array)Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00830-2 targets CXCL8/IL-8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CXCL8/IL-8 Recombinant antibody, PBS Only (Capture) 83903-2-PBS (240889F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: CXCL8/IL-8 Recombinant antibody, PBS Only (Detector) 83903-1-PBS (240889C2). 100 µg. Concentration 1 mg/mL.

Alternative CXCL8/IL-8 matched antibody pairs: MP00830-1, MP00830-3, MP50089-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00830-2, CXCL8/IL-8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83903-2-PBS. Detection antibody: 83903-1-PBS. Standard: Eg0152. Range: 0.625-40 ng/mL