

# NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00833-1

## Capture Antibody Information

**Catalog Number:**  
83882-4-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240809G11

**Reactivity:**  
human

**GenBank:**  
BC052582

**Conjugate:**  
Unconjugated

**Full name:**  
natural cytotoxicity triggering receptor 3

**Gene ID:**  
259197

## Detection Antibody Information

**Catalog Number:**  
83882-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240809C7

**Reactivity:**  
human

**GenBank:**  
BC052582

**Conjugate:**  
Unconjugated

**Full name:**  
natural cytotoxicity triggering receptor 3

**Gene ID:**  
259197

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00833-1 targets NCR3/NKp30 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NCR3/NKp30 Recombinant antibody, PBS Only (Capture) 83882-4-PBS (240809G11). 100 µg. Concentration 1 mg/mL.

Detection antibody: NCR3/NKp30 Recombinant antibody, PBS Only (Detector) 83882-2-PBS (240809C7). 100 µg. Concentration 1 mg/mL.

Alternative NCR3/NKp30 matched antibody pairs: MP00833-2, MP00833-3, MP00833-4, MP50119-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00833-1, NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83882-4-PBS. Detection antibody: 83882-2-PBS. Standard: Eg31663. Range: 0.313-40 ng/mL.