

# Mouse NCAM1/CD56 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00841-1**

## Capture Antibody Information

**Catalog Number:**  
83883-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240936B9  
**Reactivity:**  
mouse, rat  
**GenBank:**

**Conjugate:**  
Unconjugated  
**Full name:**  
neural cell adhesion molecule 1  
**Gene ID:**  
17967

## Detection Antibody Information

**Catalog Number:**  
83883-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240936G9  
**Reactivity:**  
mouse  
**GenBank:**

**Conjugate:**  
Unconjugated  
**Full name:**  
neural cell adhesion molecule 1  
**Gene ID:**  
17967

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.469-60 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00841-1 targets NCAM1/CD56 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse NCAM1/CD56 Recombinant antibody, PBS Only (Capture) 83883-2-PBS (240936B9). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse NCAM1/CD56 Recombinant antibody, PBS Only (Detector) 83883-4-PBS (240936G9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00841-1, MOUSE NCAM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83883-2-PBS. Detection antibody: 83883-4-PBS. Standard: Eg0652. Range: 0.469-60 ng/mL.