

UBAP2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00850-1

Capture Antibody Information

Catalog Number: 83885-4-PBS	Clone ID: 240890H5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ubiquitin associated protein 2
Isotype: IgG	GenBank: NM_018449	Gene ID: 55833
Purification Method: Protein A purification	Immunogen Catalog Number: Ag30572	

Detection Antibody Information

Catalog Number: 83885-3-PBS	Clone ID: 240890G11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ubiquitin associated protein 2
Isotype: IgG	GenBank: NM_018449	Gene ID: 55833
Purification Method: Protein A purification	Immunogen Catalog Number: Ag30572	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00850-1 targets UBAP2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: UBAP2 Recombinant antibody, PBS Only (Capture/Detector) 83885-4-PBS (240890H5). 100 µg. Concentration 1 mg/mL.

Detection antibody: UBAP2 Recombinant antibody, PBS Only (Detector) 83885-3-PBS (240890G11). 100 µg. Concentration 1 mg/mL.

Alternative UBAP2 matched antibody pairs: MP00850-2, MP00850-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00850-1, UBAP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83885-4-PBS. Detection antibody: 83885-3-PBS. Standard: Ag30572. Range: 0.625-80 ng/mL