

CHRNA4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00853-1

Capture Antibody Information

Catalog Number: 83915-3-PBS	Clone ID: 240986D10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: cholinergic receptor, nicotinic, gamma
Isotype: IgG	GenBank: BC111802	Gene ID: 1146
Purification Method: Protein A purification	Immunogen Catalog Number: Ag16489	

Detection Antibody Information

Catalog Number: 83915-2-PBS	Clone ID: 240986B6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: cholinergic receptor, nicotinic, gamma
Isotype: IgG	GenBank: BC111802	Gene ID: 1146
Purification Method: Protein A purification	Immunogen Catalog Number: Ag16489	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00853-1 targets CHRNA4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHRNA4 Recombinant antibody, PBS Only (Capture) 83915-3-PBS (240986D10). 100 µg. Concentration 1 mg/mL.

Detection antibody: CHRNA4 Recombinant antibody, PBS Only (Detector) 83915-2-PBS (240986B6). 100 µg. Concentration 1 mg/mL.

Alternative CHRNA4 matched antibody pairs: MP00853-2, MP00853-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00853-1, CHRNG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83915-3-PBS. Detection antibody: 83915-2-PBS. Standard: Ag16489. Range: 0.625-40 ng/mL