

PCDHGB3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00862-2

Capture Antibody Information

Catalog Number:
83879-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241052D5
Reactivity:
human
GenBank:
BC150172
Immunogen Catalog Number:
Ag15569

Conjugate:
Unconjugated
Full name:
protocadherin gamma subfamily B, 3
Gene ID:
56102

Detection Antibody Information

Catalog Number:
83879-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241052D8
Reactivity:
human
GenBank:
BC150172
Immunogen Catalog Number:
Ag15569

Conjugate:
Unconjugated
Full name:
protocadherin gamma subfamily B, 3
Gene ID:
56102

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00862-2 targets PCDHGB3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCDHGB3 Recombinant antibody, PBS Only (Capture) 83879-3-PBS (241052D5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: PCDHGB3 Recombinant antibody, PBS Only (Detector) 83879-4-PBS (241052D8). 100 μ g. Concentration 1 mg/mL.

Alternative PCDHGB3 matched antibody pairs: MP00862-1, MP00862-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00862-2, PCDHGB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83879-3-PBS. Detection antibody: 83879-4-PBS. Standard: Ag15569. Range: 0.625-40 ng/mL.