

Securin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00876-1

Capture Antibody Information

Catalog Number:
83981-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241090E11

Reactivity:
human

GenBank:
BC101834

Immunogen Catalog Number:
Ag12670

Conjugate:
Unconjugated

Full name:
pituitary tumor-transforming 1

Gene ID:
9232

Detection Antibody Information

Catalog Number:
83981-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241090G10

Reactivity:
human

GenBank:
BC101834

Immunogen Catalog Number:
Ag12670

Conjugate:
Unconjugated

Full name:
pituitary tumor-transforming 1

Gene ID:
9232

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00876-1 targets Securin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Securin Recombinant antibody, PBS Only (Capture) 83981-2-PBS (241090E11). 100 µg. Concentration 1 mg/mL.

Detection antibody: Securin Recombinant antibody, PBS Only (Capture/Detector) 83981-4-PBS (241090G10). 100 µg. Concentration 1 mg/mL.

Alternative Securin matched antibody pairs: MP00876-2, MP00876-3, MP00876-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00876-1, Securin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83981-2-PBS. Detection antibody: 83981-4-PBS. Standard: Ag12670. Range: 0.78-100 ng/mL