

# PON1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00877-2**

## Capture Antibody Information

**Catalog Number:**  
83982-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
241091B7

**Reactivity:**  
human

**GenBank:**  
BC074719

**Immunogen Catalog Number:**  
Ag12934

**Conjugate:**  
Unconjugated

**Full name:**  
paraoxonase 1

**Gene ID:**  
5444

## Detection Antibody Information

**Catalog Number:**  
83982-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
241091B2

**Reactivity:**  
human

**GenBank:**  
BC074719

**Immunogen Catalog Number:**  
Ag12934

**Conjugate:**  
Unconjugated

**Full name:**  
paraoxonase 1

**Gene ID:**  
5444

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00877-2 targets PON1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PON1 Recombinant antibody, PBS Only (Capture) 83982-2-PBS (241091B7). 100 µg. Concentration 1 mg/mL.

Detection antibody: PON1 Recombinant antibody, PBS Only (Detector) 83982-1-PBS (241091B2). 100 µg. Concentration 1 mg/mL.

Alternative PON1 matched antibody pairs: MP00877-1, MP00877-3, MP00877-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00877-2, PON1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83982-2-PBS. Detection antibody: 83982-1-PBS. Standard: Ag12934. Range: 0.781-100 ng/mL.