

HSPC159 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00879-2

Capture Antibody Information

Catalog Number: 83983-1-PBS	Clone ID: 241113D2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: galectin-related protein
Isotype: IgG	GenBank: BC062691	Gene ID: 29094
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10109	

Detection Antibody Information

Catalog Number: 83983-2-PBS	Clone ID: 241113F3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: galectin-related protein
Isotype: IgG	GenBank: BC062691	Gene ID: 29094
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10109	

Applications

Tested Applications: Cytometric bead array	Range: 0.78-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

Product Information

MP00879-2 targets HSPC159 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HSPC159 Recombinant antibody, PBS Only (Capture) 83983-1-PBS (241113D2). 100 µg. Concentration 1 mg/mL.

Detection antibody: HSPC159 Recombinant antibody, PBS Only (Detector) 83983-2-PBS (241113F3). 100 µg. Concentration 1 mg/mL.

Alternative HSPC159 matched antibody pairs: MP00879-1, MP00879-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00879-2, HSPC159 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83983-1-PBS. Detection antibody: 83983-2-PBS. Standard: Ag10109. Range: 0.78-100 ng/mL.