

Histone H1.0 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00895-3**

Capture Antibody Information

Catalog Number: 83993-2-PBS	Clone ID: 241147D11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: H1 histone family, member O
Isotype: IgG	GenBank: BC000145	Gene ID: 3005
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9982	

Detection Antibody Information

Catalog Number: 83993-1-PBS	Clone ID: 241147D2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human, mouse, rat	Full name: H1 histone family, member O
Isotype: IgG	GenBank: BC000145	Gene ID: 3005
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9982	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP00895-3 targets Histone H1.0 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Histone H1.0 Recombinant antibody, PBS Only 83993-2-PBS (241147D11). 100 µg. Concentration 1 mg/mL.

Detection antibody: Histone H1.0 Recombinant antibody, PBS Only (Detector) 83993-1-PBS (241147D2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00895-3, Histone H1.0 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83993-2-PBS. Detection antibody: 83993-1-PBS. Standard: Ag9982. Range: 0.781-100 ng/mL.