

ATG4D Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00902-3**

Capture Antibody Information

Catalog Number: 83951-1-PBS	Clone ID: 241035A7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ATG4 autophagy related 4 homolog D (S. cerevisiae)
Isotype: IgG	GenBank: BC068992	Gene ID: 84971
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10177	

Detection Antibody Information

Catalog Number: 83951-4-PBS	Clone ID: 241035G8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human, rat	Full name: ATG4 autophagy related 4 homolog D (S. cerevisiae)
Isotype: IgG	GenBank: BC068992	Gene ID: 84971
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10177	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00902-3 targets ATG4D in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ATG4D Recombinant antibody, PBS Only (Capture) 83951-1-PBS (241035A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: ATG4D Recombinant antibody, PBS Only (Detector) 83951-4-PBS (241035G8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00902-3, ATG4D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83951-1-PBS. Detection antibody: 83951-4-PBS. Standard: Ag10177. Range: 0.625-40 ng/mL