

C10orf33 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00905-2

Capture Antibody Information

Catalog Number:
83954-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240984G10

Reactivity:
human

GenBank:
BC006131

Immunogen Catalog Number:
Ag14701

Conjugate:
Unconjugated

Full name:
chromosome 10 open reading frame 33

Gene ID:
84795

Detection Antibody Information

Catalog Number:
83954-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240984F9

Reactivity:
human

GenBank:
BC006131

Immunogen Catalog Number:
Ag14701

Conjugate:
Unconjugated

Full name:
chromosome 10 open reading frame 33

Gene ID:
84795

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00905-2 targets C10orf33 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PYROXD2 Recombinant antibody, PBS Only (Capture) 83954-3-PBS (240984G10). 100 µg. Concentration 1 mg/mL.

Detection antibody: PYROXD2 Recombinant antibody, PBS Only (Capture/Detector) 83954-2-PBS (240984F9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00905-2, C10orf33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83954-3-PBS. Detection antibody: 83954-2-PBS. Standard: Ag14701. Range: 0.625-80 ng/mL.