

XBP1S Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00917-3**

Capture Antibody Information

Catalog Number: 83959-2-PBS	Clone ID: 241004E5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: X-box binding protein 1
Isotype: IgG	GenBank: BC000938	Gene ID: 7494
Purification Method: Protein A purification	Immunogen Catalog Number: Ag21703	

Detection Antibody Information

Catalog Number: 83959-1-PBS	Clone ID: 241004E1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: X-box binding protein 1
Isotype: IgG	GenBank: BC000938	Gene ID: 7494
Purification Method: Protein A purification	Immunogen Catalog Number: Ag21703	

Applications

Tested Applications: Cytometric bead array	Range: 0.156-10 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00917-3 targets XBP1S in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: XBP1S-specific Recombinant antibody, PBS Only (Capture) 83959-2-PBS (241004E5). 100 µg. Concentration 1 mg/mL.

Detection antibody: XBP1S-specific Recombinant antibody, PBS Only (Detector) 83959-1-PBS (241004E1). 100 µg. Concentration 1 mg/mL.

Alternative XBP1S matched antibody pairs: MP00917-1, MP00917-2, MP00917-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00917-3, XBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83959-2-PBS. Detection antibody: 83959-1-PBS. Standard: Ag21703. Range: 0.156-10 ng/mL