

GALNT10 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00920-1**

Capture Antibody Information

Catalog Number:
83962-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241085B9
Reactivity:
human
GenBank:
BC007224
Immunogen Catalog Number:
Ag14751

Conjugate:
Unconjugated
Full name:
UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase 10 (GalNAc-T10)
Gene ID:
55568

Detection Antibody Information

Catalog Number:
83962-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241085C2
Reactivity:
human
GenBank:
BC007224
Immunogen Catalog Number:
Ag14751

Conjugate:
Unconjugated
Full name:
UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase 10 (GalNAc-T10)
Gene ID:
55568

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00920-1 targets GALNT10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GALNT10 Recombinant antibody, PBS Only (Capture) 83962-1-PBS (241085B9). 100 µg. Concentration 1 mg/mL.

Detection antibody: GALNT10 Recombinant antibody, PBS Only (Detector) 83962-2-PBS (241085C2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00920-1, GALNT10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83962-1-PBS. Detection antibody: 83962-2-PBS. Standard: Ag14751. Range: 0.625-80 ng/mL.