

UTP18 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00941-1

Capture Antibody Information

Catalog Number:
84035-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241240C1
Reactivity:
human
GenBank:
NM_016001
Immunogen Catalog Number:
Ag34662

Conjugate:
Unconjugated
Full name:
UTP18, small subunit (SSU)
processome component, homolog
(yeast)
Gene ID:
51096

Detection Antibody Information

Catalog Number:
84035-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241240G6
Reactivity:
human
GenBank:
NM_016001
Immunogen Catalog Number:
Ag34662

Conjugate:
Unconjugated
Full name:
UTP18, small subunit (SSU)
processome component, homolog
(yeast)
Gene ID:
51096

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00941-1 targets UTP18 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: UTP18 Recombinant antibody, PBS Only (Capture) 84035-1-PBS (241240C1). 100 µg. Concentration 1 mg/mL.

Detection antibody: UTP18 Recombinant antibody, PBS Only (Detector) 84035-4-PBS (241240G6). 100 µg. Concentration 1 mg/mL.

Alternative UTP18 matched antibody pairs: MP00941-2, MP00941-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00941-1, UTP18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84035-1-PBS. Detection antibody: 84035-4-PBS. Standard: Ag34662. Range: 0.781-100 ng/mL.