

EXOC1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00942-3**

Capture Antibody Information

Catalog Number:
84007-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241058B2

Reactivity:
human

GenBank:
BC020650

Immunogen Catalog Number:
Ag2303

Conjugate:
Unconjugated

Full name:
exocyst complex component 1

Gene ID:
55763

Detection Antibody Information

Catalog Number:
84007-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241058D10

Reactivity:
human

GenBank:
BC020650

Immunogen Catalog Number:
Ag2303

Conjugate:
Unconjugated

Full name:
exocyst complex component 1

Gene ID:
55763

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00942-3 targets EXOC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EXOC1 Recombinant antibody, PBS Only (Capture) 84007-2-PBS (241058B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: EXOC1 Recombinant antibody, PBS Only (Capture/Detector) 84007-3-PBS (241058D10). 100 µg. Concentration 1 mg/mL.

Alternative EXOC1 matched antibody pairs: MP00942-1, MP00942-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00942-3, EXOC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84007-2-PBS. Detection antibody: 84007-3-PBS. Standard: Ag2303. Range: 0.313-40 ng/mL