

CCDC91 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00946-2

Capture Antibody Information

Catalog Number: 84010-2-PBS	Clone ID: 241159C8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: coiled-coil domain containing 91
Isotype: IgG	GenBank: BC028682	Gene ID: 55297
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34863	

Detection Antibody Information

Catalog Number: 84010-3-PBS	Clone ID: 241159E2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: coiled-coil domain containing 91
Isotype: IgG	GenBank: BC028682	Gene ID: 55297
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34863	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00946-2 targets CCDC91 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCDC91 Recombinant antibody, PBS Only (Capture) 84010-2-PBS (241159C8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CCDC91 Recombinant antibody, PBS Only (Detector) 84010-3-PBS (241159E2). 100 µg. Concentration 1 mg/mL.

Alternative CCDC91 matched antibody pairs: MP00946-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00946-2, CCDC91 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84010-2-PBS. Detection antibody: 84010-3-PBS. Standard: Ag34863. Range: 0.313-40 ng/mL.