

NFKB2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00950-2

Capture Antibody Information

Catalog Number:
84022-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241126F6

Reactivity:
human

GenBank:
BC002844

Immunogen Catalog Number:
Ag0662

Conjugate:
Unconjugated

Full name:
nuclear factor of kappa light polypeptide gene enhancer in B-cells 2 (p49/p100)

Gene ID:
4791

Detection Antibody Information

Catalog Number:
84022-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241126F8

Reactivity:
human

GenBank:
BC002844

Immunogen Catalog Number:
Ag0662

Conjugate:
Unconjugated

Full name:
nuclear factor of kappa light polypeptide gene enhancer in B-cells 2 (p49/p100)

Gene ID:
4791

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00950-2 targets NFKB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFKB2 Recombinant antibody, PBS Only (Capture) 84022-3-PBS (241126F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: NFKB2,p100 Recombinant antibody, PBS Only (Capture/Detector) 84022-4-PBS (241126F8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00950-2, NFkB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84022-3-PBS. Detection antibody: 84022-4-PBS. Standard: Ag0662. Range: 0.781-100 ng/mL.