

CSF2RA/CD116 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00953-2**

Capture Antibody Information

Catalog Number:
84023-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241266D10

Reactivity:
human

GenBank:
BC002635

Conjugate:
Unconjugated

Full name:
colony stimulating factor 2 receptor,
alpha, low-affinity (granulocyte-
macrophage)

Gene ID:
1438

Detection Antibody Information

Catalog Number:
84023-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241266C2

Reactivity:
human

GenBank:
BC002635

Conjugate:
Unconjugated

Full name:
colony stimulating factor 2 receptor,
alpha, low-affinity (granulocyte-
macrophage)

Gene ID:
1438

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00953-2 targets CSF2RA/CD116 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CSF2RA Recombinant antibody, PBS Only (Capture) 84023-4-PBS (241266D10). 100 µg. Concentration 1 mg/mL.

Detection antibody: CSF2RA Recombinant antibody, PBS Only (Detector) 84023-3-PBS (241266C2). 100 µg. Concentration 1 mg/mL.

Alternative CSF2RA/CD116 matched antibody pairs: MP00953-1, MP00953-3, MP01319-1, MP01319-2, MP01319-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00953-2, CSF2RA/CD116 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84023-4-PBS. Detection antibody: 84023-3-PBS. Standard: RP02335. Range: 0.625-80 ng/mL.