

PRC1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00970-2

Capture Antibody Information

Catalog Number:
84041-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241216H2
Reactivity:
human
GenBank:
BC003138
Immunogen Catalog Number:
Ag8028

Conjugate:
Unconjugated
Full name:
protein regulator of cytokinesis 1
Gene ID:
9055

Detection Antibody Information

Catalog Number:
84041-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241216H5
Reactivity:
human
GenBank:
BC003138
Immunogen Catalog Number:
Ag8028

Conjugate:
Unconjugated
Full name:
protein regulator of cytokinesis 1
Gene ID:
9055

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00970-2 targets PRC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PRC1 Recombinant antibody, PBS Only (Capture) 84041-2-PBS (241216H2). 100 µg. Concentration 1 mg/mL.

Detection antibody: PRC1 Recombinant antibody, PBS Only (Detector) 84041-3-PBS (241216H5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00970-2, PRC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84041-2-PBS. Detection antibody: 84041-3-PBS. Standard:null. Range: null.

Cytometric bead array standard curve of MP00970-2, PRC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84041-2-PBS. Detection antibody: 84041-3-PBS. Standard: Ag8028. Range: 0.625-80 ng/mL