

TREM-1/CD354 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00988-2

Capture Antibody Information

Catalog Number:
84071-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241279D7

Reactivity:
human

GenBank:
NM_018643.5

Immunogen Catalog Number:
Eg1438

Conjugate:
Unconjugated

Full name:
triggering receptor expressed on myeloid cells 1

Gene ID:
54210

Detection Antibody Information

Catalog Number:
84071-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241279E1

Reactivity:
human

GenBank:
NM_018643.5

Immunogen Catalog Number:
Eg1438

Conjugate:
Unconjugated

Full name:
triggering receptor expressed on myeloid cells 1

Gene ID:
54210

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00988-2 targets TREM-1/CD354 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TREM-1/CD354 Recombinant antibody, PBS Only (Capture) 84071-3-PBS (241279D7). 100 µg. Concentration 1 mg/mL.

Detection antibody: TREM-1/CD354 Recombinant antibody, PBS Only 84071-4-PBS (241279E1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00988-2, TREM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84071-3-PBS. Detection antibody: 84071-4-PBS. Standard: Eg1438. Range: 0.156-20 ng/mL