

CD36 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00992-1

Capture Antibody Information

Catalog Number:
84078-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241262F8
Reactivity:
human
GenBank:
NM_001001548.3

Conjugate:
Unconjugated
Full name:
CD36 molecule (thrombospondin receptor)
Gene ID:
948

Detection Antibody Information

Catalog Number:
84078-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241262C1
Reactivity:
human
GenBank:
NM_001001548.3

Conjugate:
Unconjugated
Full name:
CD36 molecule (thrombospondin receptor)
Gene ID:
948

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00992-1 targets CD36 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD36 Recombinant antibody, PBS Only (Capture/Detector) 84078-3-PBS (241262F8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD36 Recombinant antibody, PBS Only (Detector) 84078-1-PBS (241262C1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00992-1, CD36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84078-3-PBS. Detection antibody: 84078-1-PBS. Standard: Eg1070. Range: 1.563-100 ng/mL.