

TRP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00997-1

Capture Antibody Information

Catalog Number: 84082-1-PBS	Clone ID: 241242B3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: tyrosinase-related protein 1
Isotype: IgG	GenBank: BC052608	Gene ID: 7306
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10166	

Detection Antibody Information

Catalog Number: 84082-2-PBS	Clone ID: 241242D7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: tyrosinase-related protein 1
Isotype: IgG	GenBank: BC052608	Gene ID: 7306
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10166	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00997-1 targets TRP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRP1 Recombinant antibody, PBS Only (Capture) 84082-1-PBS (241242B3). 100 µg. Concentration 1 mg/mL.

Detection antibody: TRP1 Recombinant antibody, PBS Only (Detector) 84082-2-PBS (241242D7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00997-1, TRP1/TYRP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84082-1-PBS. Detection antibody: 84082-2-PBS. Standard: Ag10166. Range: 0.625-80 ng/mL.