

MYO6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01001-1**

Capture Antibody Information

Catalog Number: 84085-3-PBS	Clone ID: 241251G11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: myosin VI
Isotype: IgG	GenBank: BC146764	Gene ID: 4646
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24906	

Detection Antibody Information

Catalog Number: 84085-1-PBS	Clone ID: 241251E5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: myosin VI
Isotype: IgG	GenBank: BC146764	Gene ID: 4646
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24906	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/ml (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01001-1 targets MYO6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MYO6 Recombinant antibody, PBS Only (Capture/Detector) 84085-3-PBS (241251G11). 100 µg. Concentration 1 mg/ml.

Detection antibody: MYO6 Recombinant antibody, PBS Only (Detector) 84085-1-PBS (241251E5). 100 µg. Concentration 1 mg/ml.

Alternative MYO6 matched antibody pairs: MP01001-2, MP01001-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01001-1, MYO6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84085-3-PBS. Detection antibody: 84085-1-PBS. Standard: Ag24906. Range: 0.781-100 ng/ml.