

TCF1/TCF7 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01006-1**

Capture Antibody Information

Catalog Number:
84089-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241154B5

Reactivity:
human

GenBank:
BC048769

Immunogen Catalog Number:
Ag5792

Conjugate:
Unconjugated

Full name:
transcription factor 7 (T-cell specific, HMG-box)

Gene ID:
6932

Detection Antibody Information

Catalog Number:
84089-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241154C6

Reactivity:
human

GenBank:
BC048769

Immunogen Catalog Number:
Ag5792

Conjugate:
Unconjugated

Full name:
transcription factor 7 (T-cell specific, HMG-box)

Gene ID:
6932

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01006-1 targets TCF1/TCF7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TCF1/TCF7 Recombinant antibody, PBS Only 84089-1-PBS (241154B5). 100 µg. Concentration 1 mg/mL.

Detection antibody: TCF1/TCF7 Recombinant antibody, PBS Only 84089-3-PBS (241154C6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01006-1, TCF1/TCF7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84089-1-PBS. Detection antibody: 84089-3-PBS. Standard: Ag5792. Range: 1.56-20 ng/mL.