

RNASEH2B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01015-1**

Capture Antibody Information

Catalog Number: 84097-1-PBS	Clone ID: 241131C10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ribonuclease H2, subunit B
Isotype: IgG	GenBank: BC036744	Gene ID: 79621
Purification Method: Protein A purification	Immunogen Catalog Number: Ag33669	

Detection Antibody Information

Catalog Number: 84097-5-PBS	Clone ID: 241131H5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ribonuclease H2, subunit B
Isotype: IgG	GenBank: BC036744	Gene ID: 79621
Purification Method: Protein A purification	Immunogen Catalog Number: Ag33669	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP01015-1 targets RNASEH2B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RNASEH2B Recombinant antibody, PBS Only (Capture) 84097-1-PBS (241131C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: RNASEH2B Recombinant antibody, PBS Only (Detector) 84097-5-PBS (241131H5). 100 µg. Concentration 1 mg/mL.

Alternative RNASEH2B matched antibody pairs: MP01015-2, MP01015-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01015-1, RNASEH2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84097-1-PBS. Detection antibody: 84097-5-PBS. Standard: Ag33669. Range: 0.781-100 ng/mL.