

5HT2A Receptor Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01029-2**

Capture Antibody Information

Catalog Number: 84126-1-PBS	Clone ID: 241440C12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: 5-hydroxytryptamine (serotonin) receptor 2A
Isotype: IgG	GenBank: BC074848	Gene ID: 3356
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24056	

Detection Antibody Information

Catalog Number: 84126-3-PBS	Clone ID: 241440E10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: 5-hydroxytryptamine (serotonin) receptor 2A
Isotype: IgG	GenBank: BC074848	Gene ID: 3356
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24056	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01029-2 targets 5HT2A Receptor in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HTR2A Recombinant antibody, PBS Only (Capture/Detector) 84126-1-PBS (241440C12). 100 µg. Concentration 1 mg/mL.

Detection antibody: HTR2A Recombinant antibody, PBS Only (Detector) 84126-3-PBS (241440E10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01029-2, 5HT2A Receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84126-1-PBS. Detection antibody: 84126-3-PBS. Standard: Ag24056. Range: 0.313-40 ng/mL.