

# TSEN34 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01036-1**

## Capture Antibody Information

**Catalog Number:**  
84102-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241155F6  
**Reactivity:**  
human  
**GenBank:**  
BC020805  
**Immunogen Catalog Number:**  
Ag35120

**Conjugate:**  
Unconjugated  
**Full name:**  
tRNA splicing endonuclease 34 homolog (S. cerevisiae)  
**Gene ID:**  
79042

## Detection Antibody Information

**Catalog Number:**  
84102-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241155D12  
**Reactivity:**  
human  
**GenBank:**  
BC020805  
**Immunogen Catalog Number:**  
Ag35120

**Conjugate:**  
Unconjugated  
**Full name:**  
tRNA splicing endonuclease 34 homolog (S. cerevisiae)  
**Gene ID:**  
79042

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01036-1 targets TSEN34 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TSEN34 Recombinant antibody, PBS Only (Capture) 84102-3-PBS (241155F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: TSEN34 Recombinant antibody, PBS Only (Capture/Detector) 84102-2-PBS (241155D12). 100 µg. Concentration 1 mg/mL.

Alternative TSEN34 matched antibody pairs: MP01036-2, MP01036-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01036-1, TSEN34 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84102-3-PBS. Detection antibody: 84102-2-PBS. Standard: Ag35120. Range: 0.313-40 ng/mL