

PCSK9 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01071-1**

Capture Antibody Information

Catalog Number:
84172-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241395B8
Reactivity:
human
GenBank:
NM_174936.4

Conjugate:
Unconjugated
Full name:
proprotein convertase subtilisin/kexin type 9
Gene ID:
255738

Detection Antibody Information

Catalog Number:
84172-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241395G8
Reactivity:
human
GenBank:
NM_174936.4

Conjugate:
Unconjugated
Full name:
proprotein convertase subtilisin/kexin type 9
Gene ID:
255738

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01071-1 targets PCSK9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCSK9 Recombinant antibody, PBS Only (Capture) 84172-1-PBS (241395B8). 100 µg. Concentration 1 mg/mL.

Detection antibody: PCSK9 Recombinant antibody, PBS Only (Detector) 84172-5-PBS (241395G8). 100 µg. Concentration 1 mg/mL.

Alternative PCSK9 matched antibody pairs: MP00059-1, MP00059-2, MP00059-3, MP01071-2, MP01071-3, MP01071-4, MP50055-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01071-1, PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84172-1-PBS. Detection antibody: 84172-5-PBS. Standard: Eg0362. Range: 0.781-50 ng/mL