

PCSK9 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01071-2**

Capture Antibody Information

Catalog Number: 84172-1-PBS	Clone ID: 241395B8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: proprotein convertase subtilisin/kexin type 9
Isotype: IgG	GenBank: NM_174936.4	Gene ID: 255738
Purification Method: Protein A purification	Immunogen Catalog Number: Eg0362	

Detection Antibody Information

Catalog Number: 84172-2-PBS	Clone ID: 241395C5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: proprotein convertase subtilisin/kexin type 9
Isotype: IgG	GenBank: NM_174936.4	Gene ID: 255738
Purification Method: Protein A purification	Immunogen Catalog Number: Eg0362	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01071-2 targets PCSK9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCSK9 Recombinant antibody, PBS Only 84172-1-PBS (241395B8). 100 µg. Concentration 1 mg/mL.

Detection antibody: PCSK9 Recombinant antibody, PBS Only 84172-2-PBS (241395C5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01071-2, PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84172-1-PBS. Detection antibody: 84172-2-PBS. Standard: Eg0362. Range: 0.781-50 ng/mL